

CatchGene® Catch-cfDNA Serum/Plasma M4000 Kit

Kit Content

	2rxn	30rxn	
LV Module (with 50ml tube)	2	30	set
M Column	2	30	pcs
Collection Tube (2 ml)	4	60	pcs
Buffer AE	1.5	10	ml
Proteinase K	11	130	mg
Buffer DCL	9.6	144	ml
Buffer CW1 (concentrated)	7.5	120	ml
Buffer CW2 (concentrated)	1.68	25	ml
Elution Buffer	0.48	7.2	ml

Cat. NO.	
MC10002_M4000	2
MC10030 M4000	30

Important Notice ! "M Column" should be stored at 4°C upon arrival for long term storage. "Proteinase K" should be stored at -20°C upon arrival for long term storage.

Kit Preparation

1. Prepare 20 mg/ml Proteinase K

For 11 mg Proteinase K, please add 0.55 ml Buffer AE into tube and vortex thoroughly for dissolving For 130 mg Proteinase K, please add 6.5 ml Buffer AE into tube and vortex thoroughly for dissolving After dissolving into solvent, please store in 4°C for 6 month or -20°C for 1 year.

2. Prepare Buffer CW1

Add equal volume of 100% EtOH into Buffer CW1 (concentrated) to get Buffer CW1. After adding 100% EtOH, please check the sticker on the bottle and close the cap tightly.

3. Prepare Buffer CW2

Add equal volume of 100% EtOH into Buffer CW2 (concentrated) to get Buffer CW2. After adding 100% EtOH, please check the sticker on the bottle and close the cap tightly.

Sample Pretreatment

The half life of cfDNA in whole blood is very short. So, after sampling by blood tube, please must perform pretreatment as soon as possible.

For Whole Blood Sample

- 1. Centrifuge whole blood at 3,000 x g for 10 minute at room temperature.
- 2. Transfer upper layer to the 1.5 ml micro-centrifuge tube (not provided). Please avoid aspirating any cell debris or WBC, and intermediate layer, otherwise will co-extract gDNA form intact cell
- 3. Centrifuge at 11,000 x g for 10 min and transfer the supernatant for following extraction.

Please keep samples into -20°C or -70°C if extraction won't be performed immidiately after pretreatment.

General Protocol

- For **1ml**, **2ml** or **4 ml** Serum/Plasma sample
- 1. Add <u>50 µl (▲1ml)</u> or <u>100 µl (●2ml)</u> or <u>200 µl (■4 ml)</u> Proteinase K (20 mg/ml) into the bottom of 50 ml tube.
- Transfer ▲1ml, ●2ml or ■4 ml of sample (already centrifuged with high speed) to 50 ml tube, vortex 10 sec to let Proteinase K and samples mix thoroughly. Important! Avoid mixing Proteinase K with Buffer DCL directly. It will cause malfunction of Proteinase K and shortage of cfDNA. Please must mix Proteinase K with sample first, then add Buffer DCL.
- 4. Add <u>1 ml (▲1ml) or 2 ml (●2ml) or 4 ml (■4 ml)</u> Buffer DCL to 50 ml tube, vortex 30 sec.
- 5. 56°C incubate for 30 min, then cool down to room temperature (25°C)
- 6. Add <u>1 ml (▲1ml) or 2 ml (●2ml) or 4 ml (■4 ml)</u> 100% EtOH, vortex 15 sec.
- 7. Connect LV Module with M Column to become LV Column Module. Please refer to the illustration in next page.
- 8. Transfer all lysate into LV Column Module, centrifuge at 2,700 x g for 7 min, discard the flow-through. (If there is any remaining lysate above membrane, please centrifuge additional 2 min to let all lysate pass through the membrane.)
- 9. Add <u>3 ml (▲1ml)</u> or <u>7 ml (●2ml or ■4 ml)</u> CW1 Buffer into LV Column Module, centrifuge at 2,700 x g for 4 min, discard the flow-through. (If there is any remaining wash buffer above the membrane, please centrifuge additional 2 min to let all buffer pass through the membrane.)
- 10. Take LV Column Module out of 50 ml tube. Disconnect the M Column from the LV Module, then place the M Column on a 2 ml Collection Tube. Please refer to the illustration in next page.
- 11. Add 700 µl CW2 Buffer into spin column, centrifuge at 11,000 x g for 1 min, discard the flow-through.
- 12. Add 700 μl 100% EtOH into spin column, centrifuge at 11,000 x g for 1 min, discard the flow-through.
- 13. Add 700 µl 100% EtOH into spin column, centrifuge at 11,000 x g for 1 min, discard the flow-through.
- 14. Place spin column on a new 2 ml Collection Tube, centrifuge at 11,000 x g for 3 min to eliminate any remaining EtOH.
- 15. Place spin column on a new 1.5 ml micro-centrifuge tube. Add 10-20 μl Elution Buffer, incubation at room temperature for 5 min, and then centrifuge at 11,000 x g for 1 min for elution.

FOR RESEARCH USE ONLY



Kit Storage

- 1. M Column is able to ship at ambient temperature for 2 weeks. After receiving the kit, if you won't use them immediately, please take them out and store at 4°C for long term storage. (Do not freeze into -20 °C)
- 2. Proteinase K is lyophilized powder able to ship at ambient temperature for 2 weeks. After receiving the kit, if you won't use them immediately, please take them out and store at -20 °C for long term storage.
- Buffer, solvent and consumables, please store at 15-25 °C.
 If a precipitate has formed in Buffer DCL, please dissolve by incubating at 60°C.



Connect LV Module with M Column



Disconnect M Column from LV Column Module



FOR RESEARCH USE ONLY