

| Cat. No. | Rxn |
|----------|-----|
| MC20004 | 4 |
| MC20050 | 50 |
| MC20250 | 250 |

Kit Description

The *Catch-miRNA* Serum/Plasma kit enables purification of 19-24 nucleotides miRNA, small RNA and less than 1000 nucleotides RNA from serum/plasma or urine samples. Based on optimized reagent buffer and silica membrane column, *Catch-miRNA* Serum/plasma kit is able to get high quality and purity of miRNA, which can be used in wide range of downstream application such as qPCR, Microarray and NGS. It provides a convenient and eco-friendly protocol without using phenol or chloroform for RNA purification.

Kit Content

| | 4rxn | 50rxn | 250rxn | |
|------------------------------|------|-------|--------|-----|
| Spin Columns (with 2ml Tube) | 4 | 50 | 250 | pcs |
| Collection Tubes (2 ml) | 8 | 100 | 500 | pcs |
| Buffer RCL1 | 0.36 | 4.5 | 22.5 | ml |
| Buffer RCL2 | 0.12 | 1.5 | 7.5 | ml |
| Buffer CRW1 (concentrated) | 0.48 | 6 | 30 | ml |
| Buffer CRW2 (concentrated) | 0.96 | 12 | 60 | ml |
| RNase-Free H ₂ O | 0.96 | 12 | 60 | ml |

Kit Preparation

- 1. Prepare Buffer CRW1**
Add 4 volume of 100% EtOH into concentrated Buffer CRW1 to get Buffer CRW1.
After adding 100% EtOH, please tick the sticker on the bottle and close the cap tightly.
- 2. Prepare Buffer CRW2**
Add 4 volume of 100% EtOH into concentrated Buffer CRW2 to get Buffer CRW2.
After adding 100% EtOH, please tick the sticker on the bottle and close the cap tightly.

General Protocol

- Pipette 250 µl serum/plasma sample into 1.5 ml micro-centrifuge tube and add 75 µl Buffer RCL1. Pulse-vortexing for 10 sec, brief spin down then incubate at 25°C (room temperature) for 3 min.
- Add 25 µl Buffer RCL2, pulse-vortexing for 10 sec, brief spin down then incubate at 25°C (room temperature) for 1 min.
- Centrifuge at 11,000 x g for 3 min.
- Transfer clear supernatant to a new 1.5 ml micro-centrifuge tube, add 330 µl Isopropanol, pulse-vortexing for 10 sec then briefly spin down.
- Transfer all mixture to Spin Column (with 2ml Tube), incubate at 25°C (room temperature) for 2 min.
- Centrifuge at 11,000 x g for 1 min.
- Change a new collection tube, add 500 µl Buffer CRW1 into spin column, centrifuge at 11,000 x g for 1 min.
- Discard the flow-through, add 500 µl Buffer CRW2 into spin column, centrifuge at 11,000 x g for 1 min.
- Repeat step 8.
- Change a new collection tube, centrifuge at 11,000 x g for 3 min.
- Place the spin column into 1.5 ml Elution Tube, add 20-30 µl RNase-Free H₂O and incubate at 25°C (room temperature) for 2 min.
- Centrifuge at 11,000 x g for 1 min for elution.

FOR RESEARCH USE ONLY