

Kit Content

	2rxn	30rxn	
LV Module (with 50ml tube)	2	30	set
Spin Column	2	30	pcs
Collection Tubes (2 ml)	4	60	pcs
Buffer RC	30	185x3	ml
Buffer RL	3	45	ml
Buffer RB	3	45	ml
Buffer RW1 (concentrated)	4	58	ml
Buffer RW2 (concentrated)	0.3	5	ml
RNase-Free H ₂ O	0.5	5	ml

Kit Storage


Upon arrival,
 1. Please store **Buffer RC** at 4 °C for long term storage.
 Other buffer, solvent and consumables, please store at 15-25 °C.

Kit Preparation





- 1. Prepare Buffer RB**
 Add 2.4 volume of 100% EtOH into Buffer RB and vortex thoroughly.
 After adding 100% EtOH, please tick the sticker on the bottle and close the cap tightly.
- 2. Prepare Buffer RW1**
 Add equal volume of 100% EtOH into Buffer RW1 (concentrated) to get Buffer RW1.
 After adding 100% EtOH, please check the sticker on the bottle and close the cap tightly.
- 3. Prepare Buffer RW2**
 Add 4 volume of 100% EtOH into Buffer RW2 (concentrated) to get Buffer RW2.
 After adding 100% EtOH, please check the sticker on the bottle and close the cap tightly.

General Protocol





- Pipette 2 ml fresh whole blood sample into 50 ml tube and add 10 ml Buffer RC, mix well by inversion.
- Incubate on ice for 15 min. (Mix 2 times by inversion during incubation.)
- Centrifuge at 400 x g for 10 min at 4 °C to form a cell pellet and discard the supernatant completely.
- Add 4 ml of RC Buffer to res-suspend the cell pellet and mix well by inversion.
- Centrifuge at 400 x g for 10 min at 4 °C to form a cell pellet and discard the supernatant completely.
- Add 1400 µl Buffer RL (add 1% β-mercaptoethanol freshly), vortex vigorously for 30 sec, brief spin down then incubate at 25 °C (room temperature) for 5 min.
- Add 1400 µl Buffer RB, vortex 10 sec, brief spin down
- Connect LV Module with the Spin Column to become LV Column Module. Please refer to the illustration in next page.
- Transfer all lysate into LV Column Module, centrifuge at 2,700 x g for 2 min, discard the flow-through.
- Add 3ml RW1 Buffer into LV Column Module, centrifuge at 2,700 x g for 2 min, discard the flow-through.
- Take LV Column Module out of 50 ml tube. Disconnect the Spin Column from the LV Module, then place the Spin Column on a 2 ml Collection Tube. Please refer to the illustration in next page.
- (Optional) On column digest of DNA with DNase I (not provided).
- Add 700 µl Buffer RW1 into spin column, centrifuge at 11,000 x g for 1 min, discard the flow-through.
- Add 700 µl Buffer RW2 into spin column, centrifuge at 11,000 x g for 1 min, discard the flow-through.
- Add 700 µl 100% EtOH into spin column, centrifuge at 11,000 x g for 1 min, discard the flow-through.
- Place spin column on a new 2 ml Collection Tube, centrifuge at 11,000 x g for 3 min to eliminate any remaining EtOH.
- Place spin column on a new 1.5 ml micro-centrifuge tube. Add 30-100 µl RNase-Free H₂O, incubation at room temperature for 5 min, and then centrifuge at 11,000 x g for 1 min for elution.



Connect LV Module with the Spin Column



Disconnect Spin Column from LV Column Module



FOR RESEARCH USE ONLY